# Technovit<sup>®</sup>

# Technovit<sup>®</sup>

- D** Gebrauchsinformation
- GB** User Instructions
- F** Mode d'emploi
- I** Istruzioni per l'uso
- E** Modo de empleo

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## **GB** User Instructions Technovit 2-Component Resins

The fast and coldcuring Technovit resins in form of powder and liquid are state of the art in their field of application.

### Advice for materialographic testings:

The specimen must be clean and free of grease before embedding. For cleaning ethanole or benzine is recommended. Depending on the specimen pour a thin layer of Technovit resin into the embedding mould first. Then put in the specimen and fill up the mould until the margin. Easy shaped specimen can be placed onto the bottom of the embedding mould and doused with Technovit resin directly.

Let flow out the resin in a thin stream to allow airbubbles to pass off. Polymerized resins can be treated with all usual mechanical methods. Especial for materialographic applications the grinding and polishing products from the Heraeus Kulzer NEW LAM series achieve best results.

### Mixing

#### A homogeneous mixture of the resin is the precondition to achieve optimal properties!

Place liquid, then put in the half of the powder and mix thoroughly. Add the 2nd half of powder after and after and mix for 40 seconds. Let the dough set down for 10–15 seconds and then stir again for at least 20 seconds. Basic mixing ratio for 2-component-resins 2:1; ratio can be varied according to the requirements.

#### Technovit 5000: Place powder in a jar (20 g). Add liquid (13 ml). Do not place liquid first!

Stir thoroughly with a spatula until the powder is completely soaked with the liquid. Avoid beating movement to include airbubbles in the resin dough. A few hard taps on the worktable will help to escape air out of the resin. For mixing Technovit hardpaper cups and wooden spatulas are recommended. The heat caused by polymerization can my lead to formation of pores. To avoid excess heat do not mix more than 30 g of dough (20 g powder/10 g liquid). Carry off heat by cooling (in a vessel with cold water) is recommended. Larger specimen has to be prepared in several layers. Please be sure, that each layer is complete polymerized before applying the next one. For absolutely bubble free embeddings, use the Technomat pressure pot

### Storage advice

Keep vessels always closed. Not store above room temperature. After expiry date the materials should not be used anymore. Keep away from reach of children.

### Danger advice

- R11 Highly flammable.
- R20/21/22 Harmful by inhalation, if swallowed and in contact with skin.
- R33 Danger of cumulative effects.
- R36/37/38 Irritating to eyes, respiratory system and skin.
- R43 May cause sensitization by skin contact.
- R52/53 Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

### Safety advice

- S9 Keep container in a well-ventilated place.
- S16 Keep away from sources of ignition – No smoking.
- S23 Do not breathe gas/fumes/vapour/spray. Avoid contact with skin.
- S24 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- S33 Take precautionary measures against static discharges.
- S36/37 Wear suitable protective clothing and gloves.
- S43 In case of fire, use sand, carbon dioxide or powder extinguishing agent, never use water!
- S61 Avoid release to the environment. Refer to special instructions/ Safety data sheets.

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\* For more information click on our website [www.technovit.de](http://www.technovit.de).

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## GB

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

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