# Technovit<sup>®</sup> 7200/7210/7230 VLC

## Gebrauchsinformationen Technovit 7200/7210/7230 VLC

### Technovit 7200 VLC

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## User instructions Technovit 7200/7210/7230 VLC

### Technovit 7200 VLC

Photocuring one-component methacrylate-based resin for embedding and subsequent production of thin ground sections in medicine and dentistry.

### Technovit 7210 VLC

Photocuring one-component precision adhesive on the basis of methacrylate for precision adhesion of polymerized specimens on acrylic glass specimen holders.

### Technovit 7230 VLC

Photocuring one-component adhesive methacrylate-based, for the fixing of tissue specimens in the embedding mould.

### Properties

Technovit 7200 VLC was specially developed for the thin grinding technique. It completely penetrates hard tissue, and with this technique no decalcification of the tissue is necessary. Resins like composite tooth filling materials or bone cements are not affected by Technovit 7200 VLC.

Polymerization of the embedding medium takes place in the blue light unit. When adhering to the irradiation times indicated, temperatures will not exceed 40 °C.

Its mechanical strength makes Technovit 7200 VLC particularly suitable for the thin grinding technique, and the sawing and grinding equipment manufactured by Exakt is best suited for this method.

All customary dyeing methods may be applied with high-quality results and in a minimum of time.

### Handling

#### Fixation

Any fixation method customary in light microscopy may be used for this photocuring embedding system.

#### Dehydration

Dehydration of the tissue may be done via an ascending glycolmethacrylate or alcohol range.

#### Example

Glycolmethacrylate/Water	=	70/30	3 hours
Glycolmethacrylate/Water	=	80/20	3 hours
Glycolmethacrylate/Water	=	90/10	3 hours
Glycolmethacrylate	=	100	3 hours
Glycolmethacrylate	=	100	3 hours
Glycolmethacrylate	=	100	3 hours

#### Infiltration

Infiltration is initiated with a mixture of alcohol/Technovit 7200 VLC = 1/1 in the alcohol range or glycolmethacrylate/Technovit 7200 VLC = 1/1 in the glycolmethacrylate range, followed by infiltration with pure Technovit 7200 VLC.

If dehydration is done mechanically and the infiltration with agitation and vacuum, dehydration times of 3 hours per stage and infiltration times of 8 hours per stage will suffice for tissues of 2 to 3 mm thickness. For dehydration and infiltration without agitation and vacuum the times per stage should be multiplied by four.

Infiltration should take place in light-proof containers or equipment as contact with light will start the polymerization.

#### Embedding

For planeparallel embedding the prepared and infiltrated tissue specimen is placed in the embedding mould on a drop of fixation adhesive Technovit 7230 VLC and gently pressed down. It is important that the surface to be examined should be on the bottom of the mould. Then Technovit 7200 VLC is poured onto the specimen in such a manner that the specimen does not float and no air bubbles are entrapped.

#### Polymerization

The polymerization of the embedding takes place in the blue light unit in two stages. During the first stage the embedding medium is polymerized almost completely at a low light intensity so that the polymerization temperature does not exceed 40 °C and no cracks will form. The final polymerization of the embedding material infiltrated in the tissue is done at high light intensity during the second stage. Total polymerization time is 6 hours max.

If a maximum of 40 ml embedding medium is used at an ambient temperature of 23 °C, polymerization temperature will not exceed 40 °C.

### Precision adhesion

The precision adhesion of tissue specimens to an acrylic glass specimen holder is done with the photocuring precision adhesive Technovit 7210 VLC. The photopolymerization of the adhesive takes place in the precision adhesion press manufactured by Exakt and is finalized after 15 minutes irradiation time.

### Dyeing

Possible with all customary dyeing methods.

### Storage advice:

Keep vessels always closed. Not store above room temperature. After expiry date the materials should not be used anymore. Keep away from reach of children.

### Delivery units

Technovit 7200 VLC	1000 ml liquid
Technovit 7210 VLC	2 x 15 ml liquid
Technovit 7230 VLC	100 ml sirup

### Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets. [www.kulzer-technik.com](http://www.kulzer-technik.com)

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

\* For more information please visit [www.kulzer-technik.com](http://www.kulzer-technik.com)

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® = marchio registrato, Kulzer GmbH, Hanau
**KULZER GmbH**