Technovit® 7100

DAS SCHNEIDBARE

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User instructions Technovit 7100

System for the embedding and cutting of histological specimens

Technovit 7100* – Cold-polymerizing resin for embedding

Histoform S – Teflon embedding mould with long life

Histoform Q – Teflon embedding mould for large specimens

Histobloc – Block-holder for fixing embedding specimens

Technovit 3040 – Fast-curing resin for mounting specimens on Histobloc

Knifeholder – holder for Histoknife H, universal application in microtomes

Technovit Histoblade – Universal disposable microtome knife

Dosing devices

* cold curing resin based on hydroxyethylmethacrylate (HEMA):

Properties

Technovit 7100 is suitable for both thin and semithin sections. Its stability and elasticity allows serial sectioning with cuts of only 1 µm thickness when using a rotation microtome.

Technovit 7100 penetrates and polymerizes evenly all tissue parts of specimens for light-microscope examination in medicine, biology and material testing. (Exception: Mineralized tissue e.g. mouse or rat bones, undecalcified ⇒ v. Technovit 9100.)

If the recommended mixture ratio is observed, Technovit 7100 polymerizes at low temperatures.

Polymerization in Histoform S at room temperature (23 °C) with a maximum temperature of 37 °C, in Histoform Q with a maximum of 45 °C. Good stretching properties on water at approx. 25 °C. Staining and enzymatic reactions using standard methods are possible without removing the resin.

Processing

Specimen preparation

Fixation and dehydration of the specimens can be carried out as usually, whereby PE and glass containers can be used.

Pre-infiltration

Mixing ratio: equal parts of 96% or absolute ethanol and base liquid Technovit 7100. Specimens remain in the solution for 1–2 hours or longer, according to size.

Infiltration

1 g hardener I (= 1 bag) is solved in 100 ml base liquid (minimum 10 min.). Vacuum and agitation can be helpful. At 4 °C the preparation solution remains stable for approx. 4 weeks. The specimens are infiltrated in a sufficient amount of preparation solution for 1–12 h, depending on specimen thickness and type of tissue.

Polymerisation

1 ml hardener is added with the help of a pipette and stirred into 15 ml of preparation solution. Please use customary dosing devices e.g. Eppendorf pipette.
 1–3 ml of the solution are poured into Histoform S or Q, then the infiltrated specimens are placed in the form and positioned as required. Excess infiltration solution is removed before orientation of the sample in the form: prepare 1–3 ml, briefly swab the infiltrated samples with a fleece and orientate in the form. Time of workability (pot life) at room temperature (23 °C) is approx. 5–7 min. At room temperature (23 °C) the specimens will cure within

approx. 2 h. A thin smear layer forms on the back of the sample, which is removed with a fleece before blocking, and is in no way an impairment.

Examples

a) Embedding in Histoform S
 Preparation solution:
 100 ml Technovit 7100
 1 g hardener I (1 bag)

= Infiltration 1–12 h or longer, depending on type and thickness of tissue

Embedding:
 15 ml preparation solution
 1 ml hardener II

= curing will take approx. 2 h (consider pot life)

b) Embedding in Histoform Q
 Preparation solution:
 100 ml Technovit 7100
 1 g hardener I (1 bag)

= Infiltration 1–24 h or longer, depending on type and thickness of tissue

Embedding:
 30 ml preparation solution
 1,5 ml hardener II

= curing will take approx. 2 h (consider pot life)

Mounting

The cured specimen in the embedding mould Histoform S or Histoform Q is mounted with the help of Technovit 3040.

- Place the Histobloc in the recess of the embedding mould Histoform S or Histoform Q.
- Mix Technovit 3040 in a volume ratio of 2 parts powder to 1 part liquid to obtain a viscous liquid.
- Pour Technovit 3040 into the recess at the back of the Histobloc up to the meniscus.
- After about 10 min. the Histobloc together with the fixed specimen can be removed from the Histoform.

Cutting with Histoknife

- Loosen the clamping screws of the knifeholder
- Place the Histoknife in the knifeholder of the Kulzer microtome
- Clamp the Histoknife into the knifeholder by tightening the three screws
- Place the holder with the blade in the microtome and fasten down
- Adjust the individual cutting angle (5° is recommended)
- Cut

Further processing

The cuts are removed dry from the knife with a pair of tweezers. A cutting or stretching liquid is not necessary. After stretching the cut on water it is placed on a clean slide (no coating necessary). Drying at 60 °C for 2 h or at 37 °C longer. Staining or enzymatic reactions may now be carried out without removing the resin.

Storage advice

Keep vessels always closed. Not store above room temperature. After expiry date the materials should not be used anymore. Keep away from reach of children.

Packing sizes

Combination package

Technovit 7100
 500 ml liquid/5 x 1 g hardener I/40 ml hardener II
 Technovit 3040 yellow
 100 g powder/80 ml liquid

Individual packages

Technovit 3040 yellow 1000 g powder
 Technovit 3040 500 ml liquid

Histoblocs, 100 pcs.
 Microtom-knifeholder, 1 pcs. 17 cm
 Technovit Histoblade, 50 pcs.
 Histoform S, 1 pcs.
 Histoform Q, 1 pcs.

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets. www.kulzer-technik.com

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

* For more information please visit www.kulzer-technik.com

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® = marchio registrato, Kulzer GmbH, Hanau
KULZER GmbH