Technovit® 2021 LC FAST

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GB Instruction for Use – Technovit® 2021 LC FAST

Product description

Technovit 2021 LC FAST is an odourless, light-curing 1-component resin based on high-molecular special monomers for transparent specimen mountings. Technovit 2021 LC FAST has been developed for quick and crystal clear mounting. Free from methyl methacrylates (MMA), tetrahydrofuranes (THF) and CMR substances.

Properties

Technovit 2021 LC FAST Liquid is light-curing and is cured with non-hazardous, visible blue light in the Kulzer Technotray POWER light polymerisation unit. The unlimited pot life allows optimum positioning and infiltration of the specimen. The 1-component formulation eliminates the need for mixing and hence material losses. As Technovit 2021 LC FAST does not release gases, it does not have to be processed under an extractor hood.

Technovit 2021 LC FAST produces crystal-clear, bubble-free specimens and can be optimal processed mechanically.

Preparation of the mounting

The specimens to be mounted with Technovit 2021 LC FAST must be clean and grease-free to avoid reactions with the mounting resin. And to guarantee optimum mounting. Position the specimen on the bottom of a Kulzer mounting mould corresponding to the size of the specimen and fix it, if necessary, using Technovit 2000 LC fixing paste. Place the cylinder of the mounting mould on the floor. Fill the mould to approx. 3 mm below the rim with Technovit 2021 LC FAST. The specimen must be completely covered by the mounting resin.

Equip Technotray POWER with 6 lamps.

Mounting with infiltration

For vacuum-aided infiltration, pour Technovit 2021 LC FAST only onto the relevant area of the specimen. Then place the mounting mould together with the specimen into the vacuum set and infiltrate. After infiltration, fill the mould to approx. 3 mm below the rim with Technovit 2021 LC FAST.

Polymerisation

Place the filled mould onto the transparent drawer of the open Technotray POWER, set the timer to 5 minutes and close the drawer. The Technotray POWER starts the polymerisation process and switches off automatically at the end of the set time.

Further processing

After cooling, the specimens can be excellently ground and polished manually or automatically. The finished specimens can be easily cleaned with the usual cleaning agents. We recommend a water-based cleaning agent. Technovit 2021 is resistant to alcohol and acids, but the immersion time in alcohol-based baths should be kept as short as possible for all mounted specimens

Delivery units

Art. No. 66094513 Technovit 2021 LC FAST,

1 x 1000 ml liquid

Art. No. 66005103 Technovit 2000 LC fixing paste,

1 x 4 g

Art. No. 66060914 Technotray® POWER

light-curing unit, incl. lamp

Hazard warnings/safety advice

Please observe the instructions on the product packagings and in the safety data sheets (to be found at www.kulzer-technik.de)

Our technical advice - whether verbal, in writing or by way of trials - is given to the best of our knowledge and belief but without warranty, also with respect to any third-party property rights, and does not release you from the obligation to carry out your own tests on the products supplied by us as to their suitability for the intended processes and purposes. Application, use and processing of the products are carried out beyond our control and are therefore exclusively your responsibility. We naturally guarantee the perfect quality of our products in accordance with our general terms and conditions of sale and delivery.

Detailed information on our products can be found at www.kulzer-technik.de

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® = Registered trademark of Kulzer GmbH, Hanau
® = marca registrada, Kulzer GmbH, Hanau
KULZER GmbH